


# ORTHO Premier




*Westar Medical Products* ❖

18930 59th Avenue N.E.  
Arlington, WA 98223

[www.westarmed.com](http://www.westarmed.com)

[westar\\_medical@msn.com](mailto:westar_medical@msn.com)

T. +360 435 8900 - +800 499 3945 - Fax +425 328 1544



*Westar Medical Products* ❖

ART AND SCIENCE





# Elegance & Distinction



## Premier Orthodontic Patient Chair

Pedestal base with 5" of manual height adjustment - (shown with double articulating headrest, plush upholstery and trim)  
 Plush Upholstery 2000-060T \$3986  
 Sling Plush Upholstery 2000-060S \$3986  
 Sling Asepsis Upholstery 2000-060SA \$3492  
 Asepsis Upholstery 2000-060 \$3378  
 (Optional arms must be added with Sling upholstery)

### Standard Features:

Pedestal base with 5" manual height adjustment  
 Foot operation switch with automatic exit (return)  
 Non-articulating headrest  
 HiGlos powdercoat inside and out  
 Synchronized seat and back tilt  
 No armrests (optional add on)  
 Five year warranty  
 Quick change upholstery



## 2001 Consultation Patient Chair

Electro-mechanical - (shown with plush upholstery and trim)  
 Plush Upholstery 2000-401T \$5848  
 Sling Plush Upholstery 2000-401S \$5848  
 Sling Asepsis Upholstery 2000-401SA \$5356  
 Asepsis Upholstery 2000-401 \$5240

### Standard Features:

Electromechanical Cantilever system  
 Foot control with two programmable preset positions and automatic exit (return)  
 Left/Right access with swing up arms  
 Wheelchair accessible double articulating headrest  
 Locking swivel, 60 degree  
 HiGlos powdercoat inside and out  
 Synchronized seat and back tilt  
 Five year warranty  
 Quick change upholstery




## 2001E Economy Consultation Patient Chair

Electro-mechanical - (shown with asepsis upholstery)  
 Asepsis Upholstery 2000-017 \$4872

### Standard Features:

Electromechanical Cantilever system  
 5-Way Foot control with automatic exit (return)  
 Fixed arms  
 Wheelchair accessible double articulating headrest  
 HiGlos powdercoat inside and out  
 Synchronized seat and back tilt  
 Five year warranty  
 Quick change upholstery



Double Articulating Headrest



Ponytail Headrest



Non-articulating Headrest



Sling Backrest



Asepsis Backrest

## Orthodontic Carts and Delivery Systems

4063 Panel Mount - 3 Handpiece Automatic 4063 \$1678  
 4062 Panel Mount - 2 Handpiece Manual Wet/Dry 4062 \$ 972  
 Three way syringe with autoclavable tip  
 Standard foot control  
 Junction box accessories (no box) Air/Water shut off valves, regulators and water MSO  
 Remote mounted hanger bar (non-swivel) for panel mounts

4070 Twin Manual Control 4070 \$572  
 4071 Twin Manual Control with Vacuum 4071 \$788  
 (Junction box accessories not included with Twin units)

SP200 Ortho Cart (Casters or Fixed Base) SP200 CA or FB \$3120  
 Five drawer module, four shallow drawers and one deep drawer  
 Rear access door  
 Laminate top and sides

SP400 CPU Ortho Cart (Casters or Fixed Base) SP400 CA or FB \$3660  
 Five drawer module, four shallow drawers and one deep drawer  
 Rear access door with space for CPU  
 Keyboard pull out  
 Laminate top and sides

(upcharge may apply for non-standard colors)  
 (Cart prices do not include delivery system)

## Options

### Chair:

Double Art. Headrest (\$116 upcharge in lieu of)	7000-021	\$602
Ponytail Headrest (\$116 upcharge in lieu of)	7000-609	\$602
Non-articulating Headrest (or in lieu of)	7000-028	\$488
Moveable arms	2000-244	\$404
Stationary arms	2000-243	\$228
Upholstered arm pads	7000-096-1	\$108/set

### Delivery System:

Swivel Handpiece Holder Bar	4064	\$152
8" Handpiece Holder Bar	4003-275	\$ 54
upcharge in lieu of 4"		\$ 18
4" Handpiece Holder Bar	4003-276	\$ 36
Junction Box (no accessories)	4025-1	\$144
Water Bottle Package (Bottle only)	4003-193-1	\$186
Water Bottle System (City/Bottle switch)	4095-4	\$370
Vacuum Package (HVE, SE)	4091-7	\$212
Vacuum Package (HVE, SE, Syringe)	4091	\$344



feel free to *Explore*