



Westar Dental Equipment Price List

1-800/499-3945
Fax: 425/328-1544
westarmed.com



Westar Medical Products, Inc.

Warranty Policy

All Westar products are warranted to be free from defects in material and workmanship only. This warranty extends for a period of five (5) years unless otherwise stated. No other warranties are expressed or implied. Any and all claims against this product shall be made in writing within this time period. Remedy for correction is limited to repair or replacement of parts as determined appropriate by Westar Medical Products, Inc. This warranty shall be voided if product is misused or is not properly installed or maintained. The warranty does not cover damage resulting from the use of cleaning, disinfecting, or sterilization chemicals or processes. The warranty also does not cover o-rings, light bulbs, or other consumables, or parts sold to OEM customers. This warranty does not cover any labor or shipping charges to or from the factory. Any consequential damage will not be covered. Damage caused by accident in shipping and handling is not covered by this warranty. The Purchaser must initiate all claims against the freight carrier within twenty-four (24) hours after damaged items are received. The claim is the responsibility of the customer. Cuts and scratches in upholstery or equipment are not covered by this warranty. Upholstery packages are warranted to be free from defects in workmanship for a period of one (1) year. All lights and their parts are warranted to be free from defects in material and workmanship for a period of one (1) year. Handpiece illumination (fiber optics) will be warranted for six (6) months to be free from defect in material and workmanship only. Invoices for replacement parts may be subject to adjustment after the manufacturer has inspected the alleged defect parts. Serial numbers and purchase orders are required for warranty replacement. Most replacement parts are warranted for ninety (90) days. All Westar products are subject to continuous development; therefore, product specifications may differ from those outlined in this publication. Your Westar dealer will always have the latest information.

Order Information

Please include the following information with your purchase order:

- | | | |
|----------------------------------|----------------------------|-----------------------------|
| ❖ Catalog part number | ❖ Description of product | ❖ Desired Quantity |
| ❖ Color pattern, name and number | ❖ Shipping address | ❖ Shipment method preferred |
| ❖ Purchase Order number | ❖ Telephone and fax number | |

Terms

- ❖ All prices quoted herein are in United States dollars and are subject to change without notice.
- ❖ All prices listed in this catalog are suggested retail.
- ❖ On approved accounts terms are 1%-10 day, Net 30 from date of invoice. All other accounts will be sold on a C. O. D. or prepay basis until credit is approved.
- ❖ Shipping is F. O. B. Factory
- ❖ All orders less than \$10.00 will be charged a \$10.00 handling charge.

Returns

Merchandise may not be returned without prior authorization. Shipping must be prepaid. A restock fee of at least 15% may be charged for non-warranty related merchandise returns. Call our customer service to obtain an authorization number and additional return information. Any custom product deviating from a standard product order, including but not limited to upholstery, cannot be returned for credit consideration. Any shortages in shipment must be filed with Westar Medical Products, Inc. within five (5) working days of receipt of shipment. Claims for items damaged in transit must be filed directly with the carrier by the receiver. Authorized return items need to be returned to Westar within thirty (30) days to be eligible for credit. A copy of the original invoice or proof of purchase must accompany the authorized return. Have serial numbers and order numbers ready when requesting returns. An R. M. A. # must be written clearly on the outside of the box.

Factory Labor Charges

- ❖ Technical Labor - Retail Price: \$180 per hour Order# 1007
- ❖ Factory installation of a Westar Delivery System on a Westar chair - Retail Price: \$180 per unit Order# 1009



California Bureau of Home Furnishings and Thermal Insulation License #MFG 111108
City of Los Angeles Building and Safety Approval # M-980072

FDA Registered #9060261

All rights reserved. Printed in the U. S. A.



Table of Contents

❖ Warranty, Order Information and Terms	2
❖ Table Of Contents	3
❖ Constellation with 6000 Delivery Systems	4-5
❖ Excursion and 6000 Delivery Systems and Packages	6-7
❖ Professional Delivery Systems	8-9
❖ International with 4076 Delivery Systems	10-11
❖ 5000 Hydraulic Patient Chair and pictures of factory options	12-13
❖ Electromechanical Patient Chairs - 2001 and 2001E	14-15
❖ OS III Oral Surgery Chair	16-17
❖ OS VIII Oral Surgery Table	18-19
❖ Orthodontic Equipment, Chairs, Carts, Panel mounts	20-23
❖ Cart Delivery Systems	24-25
❖ Miscellaneous Controls	26
❖ Upholstery colors chart (printed colors are not an exact match to actual vinyl)	27
❖ WeStarlights	28
❖ Intraoral Handpiece Illumination Systems (Fiber Optics)	29
❖ Stools	30-31
❖ Upholstery Packages	32-33
❖ Junction box layout (Contoured)	34
❖ Dimensional drawings and specifications for Delivery Systems	35

Constellation Delivery Systems

Automatic control for three handpieces



Westar Assistant Constellation Package

Shown in optional Mocha Ultraleather Plush upholstery with Trim, Ultra Stools in ultraleather, and delivery mounted touchpad.


Prices as shown above: 6000 Contoured Cuspidor, shown with 2001 Patient chair, WeStarlight VIII, Ultra stools Options shown include touchpad and plush upholstery (not included in prices below) Price for package as shown above \$12,055, for touchpad add \$386

Westar Constellation Package:	Electromechanical with cuspidor	Electromechanical less cuspidor	Hydraulic with cuspidor	Hydraulic less cuspidor
Model 6000 or 6010	\$4,868	\$4,300	\$4,868	\$4,300
Model 2001 Electromechanical patient chair	\$4,762	\$4,762		
Model 5000 Hydraulic patient chair			\$5,666	\$5,666
WeStarlight VIII	\$876	\$876	\$876	\$876
D5000 Ultra Doctor Stool			\$592	\$592
A5000 Ultra Assistant Stool			\$792	\$792
D2000E Economy Doctor Stool	\$366	\$366		
A2000E Economy Assistant Stool	\$612	\$612		
Total	\$11,484	\$10,916	\$12,794	\$12,226

STANDARD FEATURES	PART NO.	MODEL	6000	6010	Euro 4035	Euro 4036
Suggested retail price						
Constellation dental unit with contoured cuspidor and assistant 's instrumentation			\$4,868		\$5,538	
Constellation dental unit with assistant's instrumentation				\$4,300		\$5,094
<u>Deluxe or Euro automatic control for three handpieces and syringe</u>			With cuspidor	Less cuspidor	With cuspidor	Less cuspidor
Standard three-way syringe, with autoclavable tip (Euro style on 4035 or 4036)	4003-120 (4003-294)		✧	✧	(✧)	(✧)
Westar arm with air brake and "S" shaped connection	4029		✧	✧	✧	✧
Stainless steel tray and anti-skid pad	4003-102		✧	✧	✧	✧
Large contoured junction box with regulators and master shut offs	5024		✧	✧	✧	✧
Standard foot control	4002-038		✧	✧	✧	✧
Purge system	4003-100		✧	✧	✧	✧
<u>Post mount utility center</u>						
Contoured porcelain cuspidor	4002-474		✧		✧	
Timed bowl rinse			✧		✧	
Momentary cup filler			✧		✧	
Telescoping arm with 2" pole clamp	4027			✧		✧
Autoclavable Premium High volume evacuator and hose	4091-1		✧	✧	✧	✧
Autoclavable Premium Saliva ejector and hose	4091-2		✧	✧	✧	✧
Self-contained water system (bottle only, see below for city/bottle option)	4003-193-1		✧	✧	✧	✧
Water outlet on utility center - 1/4" quick disconnect	4091-6		✧	✧	✧	✧
FACTORY OPTIONS			PRICE			
Additional handpiece control station	4000-050	\$168	✧	✧	✧	✧
Wet/dry foot control in lieu of standard	4003-121	\$72	✧	✧	✧	✧
Swivel instrument tray holder	4000-013-1	\$168	✧	✧	✧	✧
Touchpad for all Westar chair functions including 5 programmable presets Hydraulic chairs: Chair mounted Delivery mounted Electro-mechanical chairs: Chair mounted Delivery mounted	2000-230 2000-233 2000-231 2000-234	\$386	✧	✧	✧	✧
Upgrade to Plush upholstery with trim on chair and Ultraleather on Two stools (For details see page 32)	1006	\$554 \$334	✧	✧	✧	✧
2-way routing valve - city/bottle switch (must purchase water MSO for 6010)	4003-061	\$24	✧	✧	✧	✧
Water MSO (for use with city/bottle routing valve on 6010)	4003-141	\$84		✧		✧
Remote cup fill and bowl flush, unit-mounted	4093-1	\$124	✧		✧	
Additional High volume evacuator with hose and asepsis holder	4000-051	\$80	✧	✧	✧	✧
Additional standard syringe with tubing and asepsis holder (for assistant side)	4000-053	\$124	✧	✧	✧	✧
Additional Water/Air quick disconnect installed (1/4")	4091-6	\$34	✧	✧	✧	✧
Telescoping arm (in place of holder bar on cuspidor)	4027	\$358	✧		✧	
Internal umbilical (for Westar chairs only)	4000-054-1	\$132	✧	✧	✧	✧
Additional umbilical per foot	4000-054	\$28	✧	✧	✧	✧
Handpiece intraoral illumination system (see page 29)			✧	✧	✧	✧
Delete chair-mount bracket when adapting to other chairs	2000-226	-\$196	✧	✧	✧	✧

Excursion Delivery Systems

Rotation around the toe of Westar chairs



Westar Excursion Package

Shown in optional Silk Ultraleather Plush upholstery with Trim, Chair mounted touchpad, WeStarlight VIII, 5033 Assistant package.

Prices as shown above; Excursion with 6000 control; Shown with 2001 Patient chair with internal umbilical, Excursion WeStarlight VIII. Options shown include Excursion Assistant package #5033, touchpad, and plush upholstery (not included in prices below)

Price for package as shown above \$12,528 **without stools**.

Westar Excursion Package:	Electromechanical with 6000	Electromechanical with Euro	Hydraulic with 6000	Hydraulic with Euro
Model 5030 (Nautica) or 5037 (Euro)	\$4,392	\$4,920	\$4,392	\$4,920
Model 2001 Electromechanical patient chair	\$4,762	\$4,762		
Model 5000 Hydraulic patient chair			\$5,666	\$5,666
Excursion WeStarlight VIII (Part# 5032-245 includes pole and light) (front mount not included)	\$1,358	\$1,358	\$1,358	\$1,358
D5000 Ultra Doctor Stool	\$592	\$592	\$592	\$592
A5000 Ultra Assistant Stool	\$792	\$792	\$792	\$792
Total	\$11,896	\$12,424	\$12,800	\$13,328

STANDARD FEATURES	PART NO.	MODEL	5030	5037
Suggested retail price			\$4,392	
Excursion Delivery with 6000 Control				
Excursion Delivery with Euro Control				\$4,920
<u>Asepsis automatic control for three handpieces and syringe</u>			6000	Euro
Individual water coolant controls			✧	✧
Standard three-way syringe, with autoclavable tip - Asepsis style	4003-120		✧	
Standard three-way syringe, with autoclavable tip - Euro style	4003-294			✧
Asepsis flex arm with airlock	4000-430		✧	✧
Stainless steel tray and anti-skid pad	4003-102		✧	✧
Swivel tray and arm for Euro style delivery	4028			✧
Large contoured junction box with regulators and master shut offs	5024		✧	✧
Standard Foot Control	4002-038		✧	✧
Master on/off	4003-099		✧	✧
Self-contained water system (bottle only, see below for city/bottle option)			✧	✧
Purge system	4003-100		✧	✧
Excursion chair mounting bracket, post and hardware (Front mount)	4000-169		✧	✧
FACTORY OPTIONS		PRICE		
MODEL 5033 - Assistant's utility option: High-volume evacuator, saliva ejector, syringe,CV canister on Rear Arm	5033	\$944	✧	✧
Excursion flat screen monitor mount and arm (will need front mount if not purchased with unit) (includes 5039)	5034	\$774	✧	✧
Front mount (needed for lights or monitors purchased without unit)	4000-169	\$278	✧	✧
Excursion light post only (purchased with unit)	5032	\$482	✧	✧
Excursion light post only (not purchased with unit) includes front mount	5032-169	\$760	✧	✧
Keyboard for flat screen monitor mount	5043	\$382	✧	✧
Mouse pad tray for flat screen monitor mount	5044	\$62	✧	✧
Flat screen monitor mount with 2 inch pole clamp	5039-2IN	\$556	✧	✧
2-way routing valve - city/bottle switch (must also purchase water MSO)	4003-061	\$24	✧	✧
Water MSO (for use with city/bottle routing valve)	4003-141	\$84	✧	✧
Swivel instrument tray holder for 6000 control	4000-013-1	\$168	✧	
Touchpad for all Westar chair functions including 5 programmable presets				
Hydraulic chairs: Chair mounted Delivery mounted	2000-230 2000-233	\$386	✧	✧
Electro-mechanical chairs: Chair mounted Delivery mounted	2000-231 2000-234			
Upgrade to Plush upholstery with trim on chair and Ultraleather on Two stools (For details see page 32)	✧ 1006	\$554 \$334	✧	✧
Additional handpiece control station	4000-050	\$168	✧	
Additional High volume evacuator with hose and asepsis holder	4000-051	\$80	✧	✧
Internal umbilical (for Westar chairs only)	4000-054-1	\$132	✧	✧
Handpiece intraoral illumination system (see page 29)				
Belmont Bel-20 and Bel-50 Front Mount adapter	5040	\$556	✧	✧
Belmont Bel-20 Assistant Package	5038	\$1002	✧	✧
Belmont Bel-50 Assistant Package	5036	\$942	✧	✧

Professional Delivery Systems

Automatic control for three handpieces



Westar Cuspidor Professional Package

Shown in optional Curry Ultraleather Plush upholstery with Trim, Economy stools with ultraleather , and delivery mounted touchpad.

Prices as shown above: 6000 Professional Cuspidor, shown with 2001 Patient chair, WeStarlight VIII, Economy stools. Options shown include stainless instrument tray, 45" light pole, touchpad and plush upholstery (not included in prices below) Price for package as shown above \$10,577, for touchpad add \$386, (package also shows optional chair adapter)

Westar Professional Package:	Electromechanical with cuspidor	Electromechanical less cuspidor	Hydraulic with cuspidor	Hydraulic less cuspidor
Model 6000PC or 6000PA	\$2,700	\$2,406	\$2,700	\$2,406
Model 2001 Electromechanical patient chair	\$4,762	\$4,762		
Model 5000 Hydraulic patient chair			\$5,666	\$5,666
WeStarlight VIII	\$876	\$876	\$876	\$876
2" chair adapter, needed if ordering Westar chair	\$258	\$258	\$258	\$258
D5000 Ultra Doctor Stool			\$592	\$592
A5000 Ultra Assistant Stool			\$792	\$792
D2000E Economy Doctor Stool	\$366	\$366		
A2000E Economy Assistant Stool	\$612	\$612		
Total	\$9,574	\$9,280	\$10,883	\$10,589

9

International Delivery Systems

Automatic control for two handpieces




STANDARD FEATURES	PART NO.	MODEL	2000-020	2000-019
Suggested retail price			\$7,200	
International Package with cuspidor* (Sold as complete package only)				
International Package with telescoping arm (assistant)* (Sold as complete package only)				\$6,836
<u>Economy automatic control with two handpieces and syringe</u>			With cuspidor	Less cuspidor
Individual water coolant controls			✧	✧
Standard three-way syringe, with autoclavable tip	4003-120		✧	✧
Westar pole clamp flex arm without airlock (brake)	4030-1		✧	✧
Large contoured junction box with regulators and master shut offs	5024		✧	✧
Standard foot control	4002-038		✧	✧
Purge system	4003-100		✧	✧
<u>Post mount utility center</u>				
Water/Air quick disconnect installed (1/4") one each	4091-6		✧	
Pole mounted cuspidor with 8" porcelain bowl	4003-134		✧	
Manual bowl rinse on/off			✧	
Manual cup fill			✧	
Telescoping arm	4027			✧
Autoclavable Premium High volume evacuator and hose	4091-1		✧	✧
Autoclavable Premium Saliva ejector and hose	4091-2		✧	✧
45" Pole	2000-214		✧	✧
Westar Chair adaptor for 2" pole	2000-226		✧	✧
* Package includes: W2001E Economy Electromechanical Patient Chair with Asepsis uph., D2000E Doctor Economy Stool, WeStarlight VIII				
FACTORY OPTIONS		PRICE		
Airlock (Brake) in flex arm	4000-435	\$104	✧	✧
Additional handpiece control station	4000-050	\$168	✧	✧
Wet/dry foot control in lieu of standard	4003-121	\$72	✧	✧
Self contained water bottle	4003-193	\$168	✧	✧
City/Bottle water system	4095-4	\$210	✧	✧
Timed Bowl flush	4093	\$106	✧	✧
Additional High volume evacuator with hose and asepsis holder	4000-051	\$80	✧	✧
Additional standard syringe with tubing and asepsis holder (for assistant side)	4000-053	\$124	✧	✧
Telescoping arm (in place of holder bar on cuspidor)	4027	\$358	✧	✧
Stainless instrument tray with anti-skid pad	4003-102	\$54	✧	✧
Plastic instrument tray with anti-skid pad	4003-237	\$14	✧	✧
Additional umbilical per foot	4000-054	\$28	✧	✧
Handpiece intraoral illumination system (see page 29)			✧	✧

Note Package includes: W2001E Economy Electromechanical Patient Chair with Asepsis upholstery, D2000E Economy Doctor Stool, WeStarlight VIII, Two handpiece automatic with either cuspidor or Telescoping arm

5000 Hydraulic Patient Chair



Shown with optional Diplomat Blue Ultraleather Plush Upholstery




Estimated shipping dimensions and weights:

Hydraulic Chair -
41L x 32W x 28H - no upholstery
47L x 32W x 39H - with upholstery
330 lbs. plus upholstery 30 lbs.

Factory Options for Patient Chairs



Non-articulating headrest
7000-609 \$546
7000-028 \$442 \$104 upgrade



Double Articulating headrest
7000-021 \$546
\$104 upgrade



Florentine Burl Trim
2000-229 \$222/set



\$386

Touchpad chair mounted
Hydraulic 2000-230
Electromechanical 2000-231



Touchpad delivery mounted
Hydraulic 2000-233
Electromechanical 2000-234



Pedito upholstery has 6" shorter seat and Narrow backrest.
7005-716 \$104 upcharge

Shown in Crocus US 433 with Pony Tail headrest option, also shown with optional Twin delivery system

STANDARD FEATURES	PART NO.	MODEL	5000-001T	5000-001S	5000-001SA	5000-001
Suggested Retail Price 5000 Plush with trim hydraulic patient chair			\$6,220			
5000 Slings Plush with trim hydraulic patient chair				\$6,220		
5000 Slings asepsis hydraulic patient chair					\$5,770	
5000 Asepsis hydraulic patient chair						\$5,666
Foot Control with two programmable preset positions and two automatic exit positions (return), 7 function			✧	✧	✧	✧
Hydraulic Cantilever system			✧	✧	✧	✧
All-metal construction with HiGlos powdercoat inside and out			✧	✧	✧	✧
Safety shut-off base plate			✧	✧	✧	✧
Low voltage			✧	✧	✧	✧
Thermal overload protection			✧	✧	✧	✧
Left/Right patient access with arms that swing up			✧	✧	✧	✧
Locking swivel, 60 degree			✧	✧	✧	✧
Double-articulating headrest adjustable for height and wheelchair positioning			✧	✧	✧	✧
Synchronized seat and back tilt			✧	✧	✧	✧
Quick change upholstery (Plush style done in Ultraleather only)			✧	✧	✧	✧
Upholstered arms			✧	✧	✧	
Five year warranty			✧	✧	✧	✧
FACTORY OPTIONS	PART NO.	PRICE				
Chair mount adapter	2000-226	\$258	✧	✧	✧	✧
Arm covers for Asepsis 5000-001, upholstered	7000-096-1	\$98/set				✧
Touchpad with 5 programmable presets and two automatic exit (return) Chair mounted Delivery mounted	2000-230 2000-233	\$386	✧	✧	✧	✧
Non-articulating headrest (\$104 if purchased in lieu of articulating)	7000-028	\$442	✧	✧	✧	✧
Ponytail headrest (\$104 if purchased in lieu of articulating)	7000-609	\$546	✧	✧	✧	✧
Replacement Double articulating headrest with upholstery, (\$104 upgrade)	7000-021	\$546	✧	✧	✧	✧
Replacement upholstery, Plush with trim	7005-703	\$1,482	✧	✧	✧	✧
Replacement upholstery, Sling plush with trim	7005-706	\$1,482		✧		
Replacement upholstery, Sling asepsis	7005-709	\$1,030			✧	
Replacement upholstery, asepsis (slide on back)	7005-700	\$928				✧
Replacement headrest upholstery	7005-041	\$82	✧	✧	✧	✧
Ultraleather upgrade (not plush styling)	1005	\$416			✧	✧
Florentine Burl Trim	2000-229	\$222/set	✧	✧	✧	✧
Extended Three Year Warranty for each chair	1008	\$412	✧	✧	✧	✧

Electromechanical Patient Chairs



2001 Sling Asepsis Electromechanical Patient Chair

Shown in Capri
Blue US 514 Sling
Asepsis upholstery





2001 Plush Electromechanical Patient Chair

Shown in Sassy Ultraleather
Plush upholstery (also shown
with optional PSP Holder)

Estimated shipping dimensions and weights:

Electromechanical Chair -
41L x 32W x 28H - no upholstery
47L x 32W x 39H - with upholstery
300 lbs. plus upholstery 30 lbs.



2001E Economy Asepsis Electromechanical Patient Chair

Shown in Duration Bone asepsis upholstery



Front view of Asepsis
Backrest

STANDARD FEATURES	PART NO.	MODEL	2000-401T	2000-401S	2000-401SA	2000-401	2000-017	2000-017SA
Suggested Retail Price								
2001 Plush with trim electromechanical patient chair (available in ultraleather only)			\$5,316					
2001 Slings Plush with trim electromechanical patient chair (available in ultraleather only)				\$5,316				
2001 Slings asepsis electromechanical patient chair					\$4,868			
2001 Asepsis electromechanical patient chair						\$4,762		
2001E Economy Asepsis electromechanical patient chair							\$4,428	
2001E EconomySlings asepsis electromechanical patient chair								\$4,534
Foot Control with two programmable preset positions and automatic exit (return), 7 function			✧	✧	✧	✧		
Foot Control with automatic return, 5 function							✧	✧
Electromechanical Cantilever system			✧	✧	✧	✧	✧	✧
All-metal construction with HiGlos powdercoat inside and out			✧	✧	✧	✧	✧	✧
Thermal overload protection			✧	✧	✧	✧	✧	✧
Left/Right patient access with arms that swing up			✧	✧	✧	✧		
Stationary armrests							✧	✧
Locking swivel, 60 degree			✧	✧	✧	✧		
Double-articulating headrest adjustable for height and wheelchair positioning			✧	✧	✧	✧	✧	✧
Synchronized seat and back tilt			✧	✧	✧	✧	✧	✧
Quick change upholstery (Plush style done in Ultraleather only)			✧	✧	✧	✧	✧	✧
Upholstered arms			✧	✧	✧			
Five year warranty			✧	✧	✧	✧	✧	✧
FACTORY OPTIONS	PART NO.	PRICE						
Chair mount adapter	2000-226	\$258	✧	✧	✧	✧	✧	✧
Arm covers for Asepsis 2001 and 2001E, upholstered	7000-096-1	\$98/set				✧	✧	✧
Moveable arms (arms swing up) upgrade for economy fixed arms	2000-244	\$160					✧	✧
Touchpad with 5 programmable presets and one automatic exit (return) Chair mounted Delivery mounted	2000-231 2000-234	\$386	✧	✧	✧	✧		
Non-articulating headrest (\$104 if purchased in lieu of articulating)	7000-028	\$442	✧	✧	✧	✧	✧	✧
Ponytail headrest (\$104 if purchased in lieu of articulating)	7000-609	\$546	✧	✧	✧	✧	✧	✧
Replacement Double articulating headrest with upholstery (\$104 upgrade)	7000-021	\$546	✧	✧	✧	✧	✧	✧
Replacement upholstery, Plush with trim (available in ultraleather only)	7005-703	\$1,482	✧					
Replacement upholstery, Sling plush with trim(available in ultraleather only)	7005-706	\$1,482		✧				
Replacement upholstery, Sling asepsis	7005-709	\$1,030			✧			✧
Replacement upholstery, asepsis (slide on back)	7005-700	\$928				✧		✧
Replacement headrest upholstery	7005-041	\$82	✧	✧	✧	✧	✧	✧
Ultraleather upgrade (not plush styling)	1005	\$416	✧	✧	✧	✧	✧	✧
Florentine Burl Trim	2000-229	\$222/set	✧	✧	✧	✧	✧	✧
Extended Three Year Warranty for each chair	1008	\$412	✧	✧	✧	✧	✧	✧

OS III Oral Surgery Patient Chair



Shown with Eight way I. V. arm option


Shown in Duration Taupe upholstery

Estimated shipping dimensions and weights:

OS III Chair -
 41L x 32W x 28H - no upholstery
 47L x 32W x 39H - with upholstery
 300 lbs. plus upholstery 30 lbs.



Standard I. V. arm pillow
 Interchangeable Left and Right




Eight way articulating I. V. arm option
 Interchangeable Left and Right



Chinner headrest option
 7000-608 \$760
 upgrade \$318

STANDARD FEATURES	PART NO.	MODEL	2000-080-2
Suggested Retail Price OSIII Oral surgery electromechanical patient chair			\$6,430
Foot Control with two programmable preset positions and automatic exit (return), 7 function			✧
Electromechanical Cantilever system			✧
All-metal construction with HiGlos powdercoat inside and out			✧
Thermal overload protection			✧
Left/Right I. V. arms with one interchangeable pillow			✧
Double-articulating headrest with lock down for surgery feature adjustable for height and wheelchair positioning			✧
Synchronized seat and back tilt			✧
Quick change upholstery			✧
Upholstered arms			✧
Five year warranty			✧
FACTORY OPTIONS	PART NO.	PRICE	
Touchpad with 5 programmable presets and one automatic exit (return)			
Chair mounted	2000-231	\$386	
Delivery mounted	2000-234		
Chinner double articulating headrest (can be purchased as upgrade for \$318)	7000-608	\$760	
Eight way articulating I. V. arm	2000-220	\$604	
Leg and Body strap	2000-223	\$170	
Additional I. V. pillow only	7000-116	\$168	
Height spacer 3" (factory installed only)	2000-119	\$212	
Mobile locking base	7000-371	\$1,290	
Chair mount adapter	2000-226	\$258	
Replacement Double articulating headrest with upholstery	7000-021	\$546	
Replacement upholstery asepsis	7000-535	\$1,276	
Replacement headrest upholstery	7005-041	\$82	
Ultraleather upgrade (not plush styling)	1005	\$416	



OS VIII Oral Surgery Patient Table




Chinner headrest option
7000-608



STANDARD FEATURES	PART NO.	MODEL	2000-088
Suggested Retail Price OS VIII Oral surgery electromechanical patient chair			\$6,950
Foot Control with two programmable preset positions and automatic exit (return), 7 function			✧
Electromechanical Cantilever system			✧
All-metal construction with HiGlos powdercoat inside and out			✧
Flat when fully reclined			✧
Synchronized toe articulation			✧
Double-articulating headrest with lock down for surgery feature adjustable for height and wheelchair positioning			✧
Synchronized seat and back tilt			✧
60 degree swivel			✧
Adjustable / Removeable upholstered I. V. arms with one interchangeable pillow			✧
Five year warranty			✧
FACTORY OPTIONS	PART NO.	PRICE	
Touchpad with 5 programmable presets and one automatic exit (return) Chair mounted Delivery mounted	2000-231 2000-234	\$386	✧
Chinner double articulating headrest (can be purchased as upgrade for \$318)	7000-608	\$760	✧
Eight way articulating I. V. arm	2000-220	\$604	✧
Leg and Body strap	2000-223	\$170	✧
Additional I. V. pillow only	7000-116	\$168	✧
Chair mount adapter	2000-226	\$258	✧
Replacement Double articulating headrest with upholstery	7000-021	\$546	✧
Replacement upholstery asepsis	7005-530	\$1,276	✧
Replacement headrest upholstery	7005-041	\$82	✧
Ultraleather upgrade (not plush styling)	1005	\$416	✧

Estimated shipping dimensions and weights:

OS VIII Chair -
40L x 32W x 45H - includes upholstery
402 lbs.




Premier Orthodontic Patient Chair



Shown in optional Papyrus ultraeather plush upholstery with trim, optional double articulating headrest, and optional armrests

Estimated shipping dimensions and weights:


Premier Chair -
 41L x 32W x 28H - no upholstery
 47L x 32W x 39H - with upholstery
 240 lbs. plus upholstery 30 lbs.




STANDARD FEATURES	PART NO.	MODEL	2000-060T	2000-060S	2000-060SA	2000-060
Suggested Retail Price PremierPlush with trim Orthodontic patient chair (available in ultraleather only)			\$3,622			
PremierSlings Plush with trim Orthodontic patient chair (available in ultraleather only) (must upgrade to moveable or fixed arms)				\$3,622		
PremierSlings asepsis Orthodontic patient chair (must upgrade to moveable or fixed arms)					\$3,174	
PremierAsepsis Orthodontic patient chair						\$3,070
Foot Switches with automatic exit (return), mounted on pedestal			✧	✧	✧	✧
Electromechanical backrest			✧	✧	✧	✧
All-metal construction with HiGlos powdercoat inside and out			✧	✧	✧	✧
Pedestal base with 5" of manual height adjustment			✧	✧	✧	✧
Non-articulating headrest			✧	✧	✧	✧
Synchronized seat and back tilt			✧	✧	✧	✧
Quick change upholstery (Plush style done in Ultraleather only)			✧	✧	✧	✧
Five year warranty			✧	✧	✧	✧
FACTORY OPTIONS	PART NO.	PRICE				
Chair mount adapter	2000-226	\$258	✧	✧	✧	✧
Stationary arms	2000-243	\$206	✧	✧	✧	✧
Moveable arms (arms swing up)	2000-244	\$366	✧	✧	✧	✧
Upholstered arm covers for asepsis upholstery on added arm option	7000-096-1	\$98/set				✧
Touchpad with 5 programmable presets and two automatic exit (return) must also purchase Preset circuit board (see below)						
Chair mounted	2000-231	\$386	✧	✧	✧	✧
Preset circuit board (one preset will work with the pedestal switches, must purchase foot control or touchpad to utilize more than one preset) Also, must purchase Board Shield	7003-150-2A 7002-718	\$572 \$200	✧ ✧	✧ ✧	✧ ✧	✧ ✧
Foot Control - 7 function	7000-060-2	\$572	✧	✧	✧	✧
Non-articulating headrest (additional)	7000-028	\$442	✧	✧	✧	✧
Ponytail headrest (\$104 if purchased in lieu of standard)	7000-609	\$546	✧	✧	✧	✧
Double articulating headrest with upholstery (\$104 if purchase in lieu of standard)	7000-021	\$546	✧	✧	✧	✧
Replacement upholstery, Plush with trim (available in ultraleather only)	7005-726	\$1,482	✧			
Replacement upholstery, Sling plush with trim (available in ultraleather only)	7005-728	\$1,482		✧		
Replacement upholstery, Sling asepsis	7005-727	\$1,030			✧	
Replacement upholstery, asepsis (slide on back)	7005-725	\$928				✧
Replacement headrest upholstery	7005-041	\$82	✧	✧	✧	✧
Ultraleather upgrade (not plush styling)	1005	\$416	✧	✧	✧	✧
Florentine Burl Trim	2000-229	\$222/set	✧	✧	✧	✧
Extended Three Year Warranty for each chair	1008	\$412	✧	✧	✧	✧

Orthodontic Carts and Deliveries


ORTHO CARTS STANDARD FEATURES	PART NO.	MODEL	SP200FB	SP200CA	SP400FB	SP400CA
Standard Ortho cart with fixed base			\$2,266			
Standard Ortho cart with mobile base				\$2,266		
Ortho cart for CPU with fixed base					\$2,706	
Ortho cart for CPU with mobile base						\$2,706
Fivedrawer Module, four shallow drawers and one deep drawer			✧	✧		
Fivedrawer Module, four shallow drawers and one deep drawer					✧	✧
Rear access door			✧	✧		
Access door for CPU					✧	✧
Keyboard pullout					✧	✧
Laminate top and sides			✧	✧	✧	✧
Choice of Laminate colors, (may have upcharge)			✧	✧	✧	✧
FACTORY OPTIONS	PART NO.	PRICE				
Flat panel monitor mount (monitor not included)	FPMA103G	\$370			✧	✧
Install Panel Mount in Cart, each	1007	\$60				
PANEL MOUNT DELIVERY SYSTEMS AVAILABLE FOR ALL ORTHO CARTS	PART NO.	PRICE	4061	4062	4063	
Panel mount One handpiece manual includes one 4" hanger bar	4061	\$718				
Panel mount Two handpiece manual includes one 4" hanger bar	4062	\$882				
Panel mount Three handpiece automatic includes one 8" hanger bar	4063	\$1,524				
Three way syringe with autoclavable tip	4003-120		✧	✧	✧	
Standard foot control	4002-038		✧	✧	✧	
Junction box accessories (no box) Air/Water shut off valves, regulators and filters, Water MSO valve	4024-000		✧	✧	✧	
Remote mounted hanger bar (non -swivel)	4003-276 -4" 4003-275 -8"		✧	✧	✧	
FACTORY OPTIONS FOR PANEL MOUNT UNITS	PART NO.	PRICE	4061	4062	4063	
Vacuum Package (HVE, SE)	4091-7	\$192	✧	✧	✧	
Vacuum Package (HVE, SE, Syringe)	4091	\$312	✧	✧	✧	
Panel mount Clean water bottle (bottle water only)	4003-193-1	\$168	✧	✧	✧	
Panel mount Clean water system (city/bottle switch)	4095-4	\$212	✧	✧	✧	
Swivel remote hanger bar	4064	\$138	✧	✧	✧	
4" hanger bar	4003-276	\$32	✧	✧	✧	
8" hanger bar	4003-275	\$48	✧	✧	✧	
8" hanger bar in lieu of 4" on 4061 or 4062 upcharge	4003-275	\$16	✧	✧		
Junction box with lid (no accessories)	5025-1	\$130	✧	✧	✧	




SP200 Standard Ortho Cart



SP400 CPU Ortho Cart



4062 Two handpiece manual



4063 Three handpiece automatic



4070 Twin two handpiece manual \$516

4071 Twin with vacuum \$712



4061 One handpiece manual



4064 Swivel hanger bar option




Remote flush mounted hanger bar (non-swivel)
4" for 4061 and 4062
8" for 4063

Cart Delivery Systems

4085 Doctor's Cart - \$2852

4085 Dimensions



4089 Doctor's
and Assistant
Cart - \$3044
(add vacuum
instrumentation to
4085 Cart)


4081 Dual Purpose Cart - \$3378



4082 Assistant Cart - \$2218
(Shown with optional Stainless Tray)



4081 and 4082 Dimensions



STANDARD FEATURES	PART NO.	MODEL	4081	4082	4085	4089
Suggested Retail Price			\$3,378			
Dual purpose cart style delivery system						
Assistant's cart style delivery system				\$2,218		
Doctor's cart style delivery system					\$2,852	
Doctor's and assistant's cart style delivery system						\$3,044
Swing mount delivery for left/right convertibility	6010-1		✧			
Self-contained water system (bottle only, see below for city/bottle option)	4003-193-1		✧	✧	✧	✧
19" x 22" laminated surface - Surf White	4003-290		✧	✧		
Electrical outlet, AC	7003-390		✧	✧		
Autoclavable Premium High Volume Evacuator and hose	4003-008		✧	✧		✧
Autoclavable Premium Saliva Ejector and hose	4003-009		✧	✧		✧
Standard three-way syringe with autoclavable tip	4003-120		✧		✧	✧
Assistant's standard syringe with autoclavable tip	4000-053		✧	✧		
Contoured junction box with regulators and master shut offs	5024		✧	✧	✧	✧
Master on/off	4003-099		✧		✧	✧
Three handpiece automatic control	4003-453		✧		✧	✧
Adjustable coolant air and water for each handpiece			✧		✧	✧
Standard foot control	4002-038		✧		✧	✧
Purge system	4003-100		✧		✧	✧
Stainless steel tray with anti-skid pad	4003-102		✧		✧	✧
Water/air quick disconnect (1/4")	4091-6	\$34	✧			
FACTORY OPTIONS	PART NO.	PRICE				
Additional HVE, hose and holder	4000-051	\$80	✧	✧	✧	✧
Wet/dry foot control in lieu of standard	4003-121	\$72	✧		✧	✧
City/bottle upgrade (includes water MSO and routing valve)	4095-7	\$184	✧	✧	✧	✧
Swivel instrument tray holder	4000-013-1	\$168	✧	✧	✧	✧
Foot control rest	4000-241	\$160	✧		✧	✧
Storage organizers (3 bin)	4003-236	\$233	✧	✧	✧	✧
Storage organizers (4 bin)	4003-255	\$196	✧	✧	✧	✧
Storage bin (1 bin)	4003-256	\$28	✧	✧	✧	✧
Additional umbilical per foot (1 1/4")	4000-059	\$20	✧	✧		
Additional umbilical per foot (4085 and 4089 carts)	4000-060	\$20			✧	✧
Handpiece intraoral illumination system (see page 29)			✧		✧	✧
Cart H-frame	4080	\$508				
Cart top laminated (19" x 22") Surf White	4003-290	\$220				
Cart top laminated (14" x 14") Surf White	4003-291	\$212				
Touchpad for all Westar chair function (5 programmable presets) with 10' extention	4002-242-0					
Electro-mechanical chairs: Chair mounted	2000-230					
Delivery mounted	2000-233	\$550	✧	✧	✧	✧
Hydraulic chairs: Chair mounted	2000-231					
Delivery mounted	2000-234					

Miscellaneous Controls



4076-1 Two handpiece automatic
Control Head - \$662




4026-1 Pole Mount Cuspidor - \$874

Optional:

4091-5 Holder bar with HVE, SE, CV
Canister - \$218



4027 Telescoping Arm (shown with optional HVE
and SE with standard holders) - \$358



(Shown with optional WeStarlight VIII -
see page 28 for details.)

5050 Over the Patient
Hygiene - \$4890

5050 Standard Features:

- Two handpiece automatic with three way syringe
- Master on/off
- Purge system
- Water non-retracting coolant adjustment
- Asepsis flex arm with airlock (brake)
- Standard foot control
- Clean water system (bottle water only)
- HVE, SE, Solids collector
- Utility center
- Air/Water Quick disconnect
- 45" light pole with Chair adapter
- Contoured Junction box with accessories




4020-1 Three
handpiece au-
tomatic Control
Head - \$1082

5046 Standard Features:

- Two handpiece automatic
- Master on/off
- Purge system
- Water non-retracting coolant adjustment
- Standard foot control
- Clean water system (bottle water only)
- HVE, SE, Solids collector
- Contoured Junction box with accessories
- Rear swing style arm

5046 Rear Hygiene - \$3750





Westar Medical Products, Inc.

4470 Chennault Beach Road, Ste 2

Mukilteo, WA 98275

Ph: 425.290.3945

Toll Free: 800.499.3945

STANDARD UPHOLSTERY COLOR GUIDE*



ULTRALEATHER UPHOLSTERY COLOR GUIDE*

NOTE: WESTAR 'PLUSH UPHOLSTERY' ONLY
AVAILABLE IN ULTRALEATHER MATERIAL



* ADDITIONAL COLORS ARE AVAILABLE UPON REQUEST.



Printed colors are not an exact match for actual vinyl

WeStarlight Operator Lighting



STANDARD FEATURES	MODEL	PRICE
WeStarlight VIII Chair mounted operator light	2000-245	\$876
WeStarlight VIII Wall mounted operator light	2000-235	\$1,320
WeStarlights available in 120 volt or 220 volt (please specify)		
Dual Intensity color temperature: 4500 to 5500 degrees kelvin	✧	
Light pattern: 3.25" x 6.5"	✧	
Bulb: Quartz halogen 12v 60w	✧	
Scratch and discoloration resistant shield	✧	
Shadow free	✧	
FACTORY OPTIONS FOR WeSTARLIGHT	PART NO.	PRICE
Chair mount adapter	2000-226	\$258
45" Light pole - powdercoated light beige	2000-214	\$112
1.375 I. D. pole adapter (for WeStarlight VIII)	4002-071	\$24
1.71 I. D. pole adapter (for WeStarlight III)	7002-235	\$24
Replacement reflector	4003-358	\$452
Replacement bulb	4003-340	\$32
Replacement shield specify WeStarlight model (III or VIII)	4003-336	\$148
Replacement transformer - 120 volt 60 hz	4003-397	\$238

SYSTEM	PART NO.	PRICE	INSTALLED
Complete Non-variable One handpiece ISO-B, 5-pin System with Handpiece tubing (120v) Includes: non-variable illumination system, One ISO-B handpiece tubing (sand) Lamp module for 5-hole handpiece tubing	123215HP 1232REK 1241RSD 1242	\$434	\$518
Complete Non-variable One handpiece ISO-B, 6-pin System with Handpiece tubing (120v) Includes: non-variable illumination system, One ISO-B handpiece tubing (sand)	12321HP 1232REK 1240RSD	\$372	\$456
Features ✦ Fits any 5-hole or ISO-C (6 pin) handpiece ✦ Voltage regulated ✦ Each tubing independantly adjustable ✦ One handpiece tubing ✦ Brilliant long-life Krypton lamp ✦ Automatic activation by drive air ✦ Non-variable power pack for increased bulb life			
Complete Variable up to Three handpiece ISO-B, 5-pin System with One Handpiece tubing (120v) Includes: Variable illumination system, One ISO-B handpiece tubing (sand) Lamp module for 5-hole handpiece tubing	123235HP 1232R3V 1241RSD 1242	\$724	\$808
Complete Variable up to Three handpiece ISO-B, 6-pin System with One Handpiece tubing (120v) Includes: Variable illumination system, One ISO-B handpiece tubing (sand)	12323HP 1232R3V 1240RSD	\$662	\$746
Features ✦ Fits any 5-hole or ISO-C (6 pin) handpiece ✦ Voltage regulated ✦ Each tubing independantly adjustable ✦ Up to three handpiece tubings ✦ Brilliant long-life Krypton lamp ✦ Automatic activation by drive air ✦ Variable power pack for increased bulb life			
FACTORY OPTIONS			
Non-variable one handpiece illumination system power supply (120v) Includes: non-variable power pack and 120volt transformer	1232REK	\$206	\$252
Variable illumination system power supply (120v) Includes: variable power pack and 120 volt transformer	1232R3V	\$498	\$542
Special one handpiece control for 4069 Duo 12 O'clock, or for Power Arm style delivery Includes: Variable ISO-B illumination system One ISO-B handpiece tubing (sand) - 120" long Lamp module for 5 hole handpiece tubing	1232R2V 9SLG120 1242	\$474	\$558
ISO-B, 5-hole handpiece tubing, straight, sand (bulb not included) 5-hole tubings connect to any ISO-B 5-hole handpiece using #1242 lamp module	1241RSD	\$166	\$202
Lamp Module for all 5-hole handpiece tubings High intensity Krypton-halogen lamp, 80 hour life, mounts directly behind handpiece	1242	\$62	
ISO-C handpiece tubing, straight, sand (6-pin) ISO-C tubings connect directly to any 6-pin ISO-C handpiece. Lamp supplied within handpiece.	1240RSD	\$166	\$202
Note: Westar does not warranty system lamps. The purchase of a spare bulb is recommended.			



Art And Science

ELEGANCE & DISTINCTION Westar introduces affordable luxury. Always famous for reliability and efficiency. Westar adds beauty and comfort to the dental operator.


❖ 5000 Ultra Stools



"Elegance and Distinction - Westar introduces affordable luxury. Always famous for reliability and efficiency, Westar adds beauty and comfort to the dental operator. The classic look of contoured lines and leather will move you to the forefront of modern operator design."



❖ 2000 Economy Stools



STANDARD FEATURES	PART NO.	MODEL	D5000	A5000	D2000E	A2000E
Ultra doctor's stool with contoured seat	2000-104-1		\$592			
Ultra assistant's stool with 16" round seat	2000-105-1			\$792		
Economy doctor's stool	2000-103-1				\$366	
Economy assistant stool	2000-102-1					\$612
Ultra doctor's stool with 16" round seat	2000-104-3		\$592			
Ultra assistant's stool without body support	2000-105-3			\$624		
<i>Ultra stools are not sold without backs due to the tilt feature of the seat and back</i>						
Economy doctor's stool with 16" round seat	2000-103-2				\$366	
Economy doctor's stool without back with round seat	2000-103-3				\$304	
Economy assistant's stool without body support	2000-102-2					\$494
Economy assistant's stool with backrest	2000-102-3					\$778
Five point 21" base - powder coated cast aluminum			✧	✧		
Five point 22" molded plastic base					✧	✧
Synchronized seat and back tilt with quick adjustment paddles			✧	✧		
Vertical back adjustment			✧	✧		
Dual wheel, non-marring casters			✧	✧		
Dual wheel barrel casters					✧	✧
Five year warranty			✧	✧	✧	✧
19" wide contoured seat			✧		✧	
16" wide contoured back			✧	✧		
Horizontal and vertical back adjustment with knobs					✧	
Swing body support with adjustable height				✧		✧
Ratcheting quick release for entry and exit				✧		
Adjustable for right or left hand use				✧		✧
Adjustable foot ring				✧		✧
Height adjustment			20" to 25"	24.5" to 31"	17" to 22"	22.5" to 30"
FACTORY OPTIONS	PART NO.	PRICE				
Ultraleather upholstery upgrade	1006	\$168	✧	✧	✧	✧
"D" seat in lieu of round	8003-199	\$138		✧		✧
Assistant's cylinder (Ultra) 45 mm	8000-017-1	\$86 or in lieu of	✧			
Doctor's cylinder (Ultra) 45 mm	8000-018-1	\$86 or in lieu of		✧		
Assistant's cylinder (Economy) 50 mm	8000-017	\$86 or in lieu of			✧	
Doctor's cylinder (Economy) 50 mm	8000-018	\$86 or in lieu of				✧
Adjustable foot ring (Ultra)	8003-214	\$98	✧			
Adjustable foot ring (Economy)	8000-003	\$98			✧	
Replacement metal torso knob (for prior style swing and stand tube)	7002-342	\$24		✧		✧
Large backrest adjustment knob (specify grey or tan)	8003-202	\$14	✧	✧		
Ratchet mechanism for Deluxe and old style Ultra stools	8000-021	\$344		✧		
Casters for plastic base (specify color)- priced per caster, 5 needed per base	8000-002	\$8			✧	✧
Casters for metal base - non-marring - priced per caster, 5 needed per base	8000-033	\$10	✧	✧		
Plastic Ultra stool back cover	8002-020	\$106	✧	✧		

Upholstery Package Components for Field Installation

Pelton & Crane Chairman	\$1176	9000-001
Pelton & Crane Coachman	\$1176	9000-011
Adec 1005 (Snap)	\$1176	9000-010
Adec 1005 (Bolt)	\$1176	9000-010-1
Belmont/Healthco Celebrity (Black)	\$1284	9000-013
Belmont/Healthco Celebrity (White)	\$1284	9000-013-1
Adec 1020	\$1246	9000-019

W2000 (Screw on style Back)	\$928	7005-091
Westar2001/5000 Asepsis(Slide on Back)	\$928	7005-700
Westar 2001/5000 Plush (Slide on Back)	\$1482	7005-703
Westar 2001/5000 Sling Plush Style	\$1482	7005-706
Westar 2001/5000 Sling Asepsis Style	\$1030	7005-709
Premier Asepsis Style	\$928	7005-725
Premier Plush Style	\$1482	7005-726
Premier Sling Style	\$1030	7005-727
Premier Plush Sling Style	\$1482	7005-728
OS III	\$1276	7000-535
OS VIII	\$1276	7005-530
Pedo	\$1032	7005-716

Pelton & Crane Front Row Doctor	\$352	9000-002
Pelton & Crane Front Row Assistant	\$352	9000-003
Pelton & Crane Act II Doctor	\$352	9000-004
Pelton & Crane Act II Assistant	\$352	9000-005
Adec 1600 Doctor	\$352	9000-006
Adec 1620 Assistant	\$352	9000-007

Westar Deluxe Doctor	\$352	8000-100
Westar Deluxe Assistant	\$352	8000-101
Westar Economy Doctor	\$352	8000-103
Westar Economy Assistant	\$352	8000-102
Westar Ultra Doctor	\$352	8000-104
Westar Ultra Assistant	\$352	8000-105
Stool Recovery	\$364	9000-014

Front Row Dr/Asst. Seat	\$244	9002-016
Front Row Backrest	\$244	9002-017
Front Row Body Support	\$244	9002-018
Act II Dr/Asst Seat	\$244	9002-022
Act II Backrest	\$244	9002-023
Act II Body Support	\$244	9002-024
Adec 1600/1620 Seat	\$244	9002-027
Adec 1600 Backrest	\$244	9002-028
Adec 1620 Body Support	\$244	9002-029

Westar Ultra Doctor Seat	\$244	8005-002
Westar Ultra Back	\$244	8005-001
Westar Ultra 16" Round	\$244	8005-003
Westar Ultra Torso	\$244	8005-004
Westar Doctor Contoured Seat	\$244	8003-089
Westar Stool Back	\$244	8003-069
Westar 15" Round Seat	\$244	8003-039
Westar Torso(Body support) new style	\$244	8005-004
Westar Torso (Body support) old style	\$244	8003-119
Westar "D" shaped seat	\$244	8003-199
Child seat wedge	\$270	2000-246

Ultraleather, Luxor, or Finesse Vinyl Upcharges (Not plush styling on upholstery):

Complete Chairs	\$416	1005
Stools and Component Parts	\$168	1006

Chair Components:

Chairman/Coachman Headrest Pillow	\$202	9003-015
Chairman Telescoping Headrest	\$190	9002-007
Coachman Telescoping Headrest	\$190	9003-071
Adec1005 Headrest	\$190	9002-025
Celebrity Headrest	\$216	9003-096
Chairman/Coachman Artic. H/R (Bolt)	\$270	9003-114
Chairman/Coachman Art. H/R(Velcro)	\$270	9003-110
Westar Articulating H/R cushion	\$82	7005-041

Chairman Arm (each arm)	\$134	9002-005
Coachman Arm (each arm)	\$134	9003-076
Coachman Arm Support (each elbow)	\$134	9003-074
Chairman Backfiller (Right)	\$76	9003-007
Chairman Back filler (Left)	\$76	9003-008
Adec Arm slings (Snap back)	\$190	9002-039
Adec Arm slings (Bolt back)	\$190	9002-040

OS Arm (each arm pad, no pillow)	\$162	7000-115
OS I. V. pillow only (interchangeable L/R)	\$168	7000-116
Eight way I. V. arm upholstery	\$162	7000-117

Chairman Backrest	\$622	9000-018
Coachman Backrest	\$622	9003-069
Adec Back(Snap) (incl. slings)	\$622	9003-051
Adec Back(Bolt) (incl. slings)	\$622	9003-075
Celebrity (Black) Backrest	\$622	9003-091
Celebrity (Hybrid) Backrest	\$622	9003-091-2
Celebrity (White) Backrest	\$622	9003-111
Adec 1020 Backrest	\$622	9003-113
W2000 Backrest (Screw on style)	\$622	7005-001
Westar Back Asepsis Slide on	\$622	7005-060
Westar Back Plush	\$1136	7005-704
Westar Back Plush Sling	\$1136	7005-707
Westar Back Asepsis Sling	\$726	7005-715

Chairman Seat	\$812	9003-016
Coachman Seat	\$824	9003-070
Adec 1005 Seat	\$824	9003-052
Celebrity Seat	\$920	9003-092
Adec 1020 Seat	\$920	9003-109
Westar Seat Asepsis	\$812	7005-061
Westar Seat Plush	\$1190	7005-705
Pedo seat	\$892	7005-717
OS III Seat and old style screw on W*	\$812	7005-021

Chairman Toecover	\$76	9002-003
Coachman Toecover	\$76	9003-073
Adec Toecover	\$76	9002-026
Celebrity Toecover	\$76	9003-120
Westar Toecover	\$76	7000-065
OS III Toecover	\$76	7000-530

Misc. Vinyl per yard	\$60	5003-003
Misc. Ultraleather or Luxor per yard	\$270	5003-338

Estimated shipping dimensions and weights

Chair Upholstery -	47L x 11W x 23H	30 lbs. (70lbs U.P.S)
Stool Upholstery -	20L x 20W x 8H	11 lbs.

Pelton & Crane Front Row (A1) and Act II Stools
 Retail Price: **\$352**
 F/R Dr., Order# 9000-002; F/R Asst., #9000-003
 Act II Dr., Order# 9000-004; Act II Asst., # 9000-005



Pelton & Crane Chairman
 Retail Price: **\$1176**
 Order# 9000-001

Pelton & Crane Coachman
 Retail Price: **\$1176**
 Order# 9000-011


Adec 1020
 Retail Price: **\$1284**
 Order# 9000-019

Adec 1005
 Retail Price: **\$1176**
 Snap Back Order# 9000-010
 Bolt Back Order# 9000-010-1

Belmont/Healthco Celebrity
 Retail Price: **\$1284**
 Black Back Order# 9000-013
 Hybrid Back Order# 9000-013-2
 White Back Order# 9000-013-1

WESTAR ❖
Medical Products




Adec 1600 and 1620 Stools
 Retail Price: **\$352**
 1600 Dr., Order# 9000-006
 1620 Asst., Order# 9000-007









- NOTES:
- 1.) LOCAL REGULATIONS PROVIDE THAT LICENSED PLUMBERS AND ELECTRICIANS SHALL INSTALL UTILITIES
 - 2.) MAKE SURE ALL PLUMBING CONFORMS TO PREVALUING LOCAL CODES
 - 3.) PLUMBING CENTRAL VACUUM TO THE UTILITY CENTER SHOULD BE SPECIFIED BY THE CENTRAL VACUUM SUPPLIER AND TERMINATED IN UTILITY CENTER WITH 3/4" O.D. TUBE PERPENDICULAR TO THE FLOOR, SIMILAR TO THE DRAIN CONNECTION
 - 4.) GRAVITY DRAIN HAS 3/4" NOMINAL PIPE PROTRUDING 1" FROM FINISHED FLOOR
- NOTE: PLATE TRAP IN LINE TO CONFORM WITH LOCAL CODES, SUPPLIED BY CONTRACTOR

Dimensions and Specifications

Excursion Delivery Systems



Westar over the patient delivery systems



Estimated shipping dimensions and weights:

Hydraulic Chair -	41L x 32W x 28H	330 lbs. plus upholstery 30 lbs.
Electromechanical Chair -	41L x 32W x 28H	300 lbs. plus upholstery 30 lbs.
OS VIII	40L x 32W x 45H	402 lbs. including upholstery
Ortho Chair -	41L x 32W x 28H	300 lbs. plus upholstery 30 lbs.
Upholstery -	47L x 11W x 23H	30 lbs.
Control Head -	37L x 23W x 13H	38 lbs.
Cuspidor -	37L x 23W x 13H	49 lbs.
Assistant -	37L x 23W x 13H	36 lbs.
International	37L x 23W x 13H	52 lbs cuspidor (Asst. 57 lbs including CA)
42" Pole -	48L x 3W x 3H	5 lbs.
Chair adaptor -	17L x 7W x 7H	15 lbs.
Light -	33L x 17W x 14H	23 lbs.
Excursion	35L x 23W x 23W	90 lbs.
Excursion with Light post	47L x 33W x 18H	105 lbs.
D5000 - UPS: 23L x 23W x 19H LTL: 23L x 23W x 35H		31 lbs
A5000 - UPS: 23L x 23W x 19H LTL: 23L x 23W x 35H		40 lbs
D2000E - UPS: 23L x 23W x 19H LTL: 23L x 23W x 35H		25 lbs
A2000E - UPS: 23L x 23W x 19H LTL: 23L x 23W x 35H		31 lbs

Westar Medical Products, Inc.
Phone: 425/290-3945 ♦ Watts: 800/499-3945 ♦ Fax: 425/328-1544
4470 Chennault Beach Road, Suite 2 ♦ Mukilteo, WA 98275
Email: westar_medical@msn.com ♦ Website: www.westarmed.com