



*"Art is the expression of the profoundest thoughts in the simplest way"*

---- EINSTEIN



Phone: 425/290-3945 ♦ Watts: 800/499-3945 ♦ Fax: 425/328-1544  
4470 Chennault Beach Road, Suite 2 ♦ Mukilteo, WA 98275  
Email: [westar\\_medical@msn.com](mailto:westar_medical@msn.com) ♦ Website: [www.westarmed.com](http://www.westarmed.com)



## Premier Orthodontic Patient Chair

(Shown with D5000 Ultra Doctor Stool and plush upholstery including trim)

Premier Ortho Chair with Plush Upholstery and Trim	2000-060T	\$3622
Premier Ortho Chair with Sling Plush Upholstery and Trim	2000-060S	\$3622
Premier Ortho Chair with Sling Asepsis Upholstery	2000-060SA	\$3174
Premier Ortho Chair with Asepsis Upholstery	2000-060	\$3070

(Plush upholstery available in Ultraleather only)



Florentine Burl Trim 2000-229 \$222/set




Standard Fixed Headrest 7000-028 \$442

Ponytail Headrest 7000-609 \$546



Touchpad Delivery or Chair Mounted \$386  
(Part Numbers found on back cover)



## 2001 Electro-mechanical Patient Chair

(Shown with plush upholstery including trim)

2001 Electro-mechanical with Plush Upholstery and Trim	2000-401T	\$5316
2001 Electro-mechanical with Sling Plush Upholstery and Trim	2000-401S	\$5316
2001 Electro-mechanical with Sling Asepsis Upholstery	2000-401SA	\$4868
2001 Electro-mechanical with Asepsis Upholstery	2000-401	\$4762

(Plush upholstery available in Ultraleather only)



Sling Upholstery



**ECONOMY WITHOUT COMPROMISE** Craftsmanship of the highest order is a time honored tradition at Westar. The latest technology is combined with the skill of craftsmen to create hallmark products which are responsive, and assured to bring you affordable luxury. This makes Economy Without Compromise the tradition of choice.



2001E Economy Patient Chair  
(Shown with optional International Package)  
(Shown with Asepsis upholstery)

2001E Economy Electro-mechanical with Asepsis Upholstery  
2000-017 \$4428




(Shown with optional  
Eight-way I. V. Arm)


OS III Oral Surgery Chair with Standard Upholstery  
2000-080-2 \$6430



Chinner Headrest option



Eight Way Articulating I. V. Arm option  
(Interchangable Left and Right)



Standard I.V. Arm Pillow  
(Interchangable Left and Right)



5000 Hydraulic Patient Chair  
(Shown with plush upholstery including trim)

5000 Hydraulic with Plush Upholstery and Trim	5000-001T	\$6220
5000 Hydraulic with Sling Plush Upholstery and Trim	5000-001S	\$6220
5000 Hydraulic with Sling Asepsis Upholstery	5000-001SA	\$5770
5000 Hydraulic with Asepsis Upholstery	5000-001	\$5666
(Plush upholstery available in Ultraleather only)		



OS VIII Oral Surgery Table  
(shown with optional Chinner Headrest)

OS VIII Oral Surgery Table with Standard Upholstery  
2000-088 \$6950

### 5000 Hydraulic and 2001 Standard Features:

Wheelchair accessible double articulating headrest  
7 way foot control with two programmable preset positions and automatic return (exit)  
Metal construction with HiGlos powdercoat inside and out  
Quick Change Upholstery  
Synchronized seat and lumbar with back recline  
Thermal overload protection  
Low voltage  
Left/Right patient access with arms that swing up  
Locking 60 degree swivel  
Three year warranty  
Lifetime bushings throughout all pivot points

### Additional features included with 5000 Hydraulic Patient Chair:

Hydraulic cantilever system  
Safety shut-off plate

### Additional features included with 2001 Patient Chair:

Electro-mechanical cantilever system

### Premier Ortho Patient Chair Standard Features:

Non-articulating headrest  
Base mounted Foot switch with manual operate and automatic return (exit)  
Pedestal base with 5" manual height adjustment  
Metal construction with HiGlos powdercoat inside and out  
Quick Change Upholstery  
Synchronized seat and lumbar with back recline  
Thermal overload protection  
Low voltage  
Three year warranty  
Lifetime bushings throughout all pivot points  
No arms

### Optional features:

Fixed arms	2000-243	\$206
Movable Arms ( <i>arms swing up</i> )	2000-244	\$366
Upholstered covered arm pads	7000-096-1	\$98/set

### Optional Features for all chairs:

Ponytail headrest ( <i>\$104 upcharge if in lieu of</i> )	7000-609	\$546
Non-articulating headrest ( <i>Can be purchased in lieu of</i> )	7000-028	\$442
Double articulating headrest ( <i>\$104 upcharge if in lieu of</i> )	7000-021	\$546
Touch-pad chair control with 5 programmable presets ( <i>Upgrade for 2001E Economy chair required</i> )		

Touchpad for chair control

5000:	Chair mount	2000-230	\$386
	Delivery mount	2000-233	\$386
2001/OS III:	Chair mount	2000-231	\$386
	Delivery mount	2000-234	\$386
Chair Adapter for 2" pole		2000-226	\$258
45" Light Pole		2000-214	\$11
1.375 I. D. pole adapter		4002-071	\$24
WeStarlight VIII		2000-245	\$876
Internal umbilical		4000-054-1	\$132
Upholstered arms for non-plush upholstery		7000-096-1	\$98/set
Ultraleather upgrade - chairs ( <i>not plush styling</i> )		1005	\$416
Ultraleather upgrade - stools and small parts		1006	\$168
Florentine Burl Trim		2000-229	\$222/set
Extended Three Year Warranty		1008	\$412
Child Seat Wedge		2000-246	\$278
Pedo upholstery ( <i>Asepsis style Narrow back 6" shorter seat</i> )		7005-716	\$928

### 2001 Economy Standard Features:

Wheelchair accessible double articulating headrest  
Electro-mechanical cantilever system  
5 way foot control with auto-return  
Metal construction with HiGlos powdercoat inside and out  
Quick Change Upholstery  
Synchronized seat and lumbar with back recline  
Thermal overload protection  
Low voltage  
Fixed arms  
Three year warranty  
Lifetime bushings throughout all pivot points

### OS III Oral Surgery Chair Standard Features:

Left/Right upholstered I. V. arms with one interchangeable pillow  
Wheelchair accessible double articulating headrest  
Adjustable / Removable arms  
7 way foot control with two programmable preset positions and automatic return (exit)  
Metal construction with HiGlos powdercoat inside and out  
Quick Change Upholstery  
Synchronized seat and lumbar with back recline  
Thermal overload protection  
Low voltage  
Three year warranty  
Lifetime bushings throughout all pivot points  
Electro-mechanical cantilever system  
Floor to seat: 13" to 31"

### OS VIII Oral Surgery Table Standard Features:

Left/Right upholstered I. V. arms with one interchangeable pillow  
Wheelchair accessible double articulating headrest  
Flat when fully reclined  
Toe articulation  
Locking 60 degree swivel  
Narrow back upholstery  
Adjustable / Removable arms  
7 way foot control with two programmable preset positions and automatic return (exit)  
Metal construction with HiGlos powdercoat inside and out  
Quick Change Upholstery  
Synchronized toe with back recline  
Thermal overload protection  
Low voltage  
Three year warranty  
Lifetime bushings throughout all pivot points  
Electro-mechanical cantilever system  
Floor to seat: 20" to 35"  
Floor to Headrest: 43" to 58"  
Fully reclined (headrest to end of toe): 60"

### Optional Features OS III and OS VIII:

Chinner headrest	7000-608	\$760
In Lieu of Double Articulating		\$308
Eight way articulating I. V. arm board	2000-220	\$604
3" Height adjustment spacer	2000-119	\$212
Leg/Body Strap set	2000-223	\$170
Mobile Base ( <i>OS III only</i> )	7000-371	\$1290
Additional I. V. Arm Pillow only	7000-116	\$168